

Co-transfection of ARPE-19 Cells by Using Biontex K2® Transfection System

Jihyun Hwang and Fengyi Liang

Department of Anatomy, Yong Loo Lin School of Medicine, National University of Singapore, Singapore 117597

This experiment was carried out to test transfection efficiency of the K2 transfection kit.

Cell culture

ARPE-19 cells were cultured in 24-well plates (Thermo-Scientific) in antibiotic-free DMEM:F12 medium with 10% fetal bovine serum (Life Technologies). Cells were seeded at a density of 1×10^5 cells/well.

Transfection

Transfection was carried out when cells reached 90% confluency, normally one day after seeding.

Five μ l K2® Multiplier was added drop-wise into a well containing ARPE19 cells and 500 μ l DMEM-F12-FBS, and mixed by gently swaying the plates. No multiplier was added to the control well. Cells were incubated in the K2® Multiplier mixture (or normal medium for the control) for 2 hours.

At about 20 min before transfection, prepare the following amount of transfection mixture for each well. First, 0.6- μ g plasmid DNA and 1.2- μ l K2® Transfection reagent were each added to 30 μ l serum-free DMEM:F12 medium in separate tubes and mixed by gently pipetting 1-2 times. Our DNA plasmids in the current case were mCherry-C1 vector (0.15 μ g) (Clontech) and pXJ40-Flag-Juxt看odin (JN) (0.45 μ g). Then the DNA-medium mixture was added to the transfection reagent-medium mixture, and mixed by gently pipetting 1-2 times. This DNA-lipid complex solution was incubated at room temperature for 15 minutes before being added drop-wise to the wells of cells containing DMEM-F12-FBS (500 μ l) and K2 multiplier (0 μ l or 5 μ l). After mixing by gentle swaying, the plate was put into 37°C, 5% CO₂ incubator. The culture medium was changed 24 hrs after to DMEM:F12 medium (Life Technologies) with 10% fetal bovine serum. Transfection efficiency could be checked by viewing mCherry signal under fluorescent microscope (Nikon Eclipse TE2000-S). Cells were then cultured for an additional 2 hours.

Immunofluorescence staining

To further confirm the transfection efficiency, immunofluorescence staining was conducted. The cells were washed briefly in PBS (phosphate buffered saline, 0.01M, pH7.4), fixed in 3% paraformaldehyde for 30 min, washed again in PBS, and pre-incubated in PBS-T-NGS (0.3% Triton X-100, 6% Normal Goat Serum). They were then incubated in mouse anti-FLAG tag antibody (1:1000 in PBS-T-NGS, Sigma) solution overnight. The secondary antibody Alexa Fluor 488 Goat anti-mouse IgG (1:400 in TBS-

T-NGS, Life Technologies) was used to detect bound primary antibody. After coverslipping, the cells were viewed and photographed at 10x magnification under an inverted microscope attached with a digital camera (Nikon).

Results and Discussion

Cell survival appears to be at around 80% while a co-transfection efficiency of 15-20% was observed. This efficiency might increase considerably if transfection of a single plasmid instead of co-transfection of 2 plasmids, were carried out. The K2® Multiplier does not seem to significantly increase cell viability or transfection efficiency in ARPE-19 cells.




Figure 1. Co-transfection of ARPE-19 cells by using K2 transfection kit (but no addition of the K2 multiplier). Juxtalin-transfected cells were immuno-stained green (**A**). mCherry is red (**B**). Panel **C** shows the merge of **A** and **B**. Cell nuclei were stained by DAPI (blue) in **C**. Panel **D** shows phase-contrast image of the field. Co-transfection efficiency in this case is around 18%.