

Transfection of HepG2 cells with METAFECTENE

Andreas Kampkötter, Yvonne Chovolou, Wim Wätjen and Regine Kahl

Zentrum für Pharmakologie und Toxikologie, Institut für Toxikologie, Heinrich-Heine Universität Düsseldorf, D-40225 Düsseldorf, Germany, kampkoet@uni-duesseldorf.de

The import of DNA in cells of higher organisms (transfection) is a suitable method to investigate biological and medical questions. Under proper requirements the imported DNA is expressed and therefore the transfected cells can be used to investigate the properties of exo- and endogenous proteins, to analyse intracellular processes, to determine the activity of gene promoters and to work on many other purposes. The requirements to a transfection method are a low consumption of DNA, a low cytotoxicity and of course a high efficiency.

For the transfection of the human hepatoma cell line HepG2 a non-liposomal lipid, an activated dendrimer and METAFECTENE were compared with respect to transfection efficiency. The HepG2 cells were cultivated in RPMI medium containing 10% FCS, penicillin and streptomycin and then transfected with the plasmid pSV β -Gal (Promega) that encodes β galactosidase. The transfections were performed according to the standard protocols given in the manuals without optimization. After an incubation of 24 to 48 hours the transfection efficiency was photometrically determined with a β -Gal assay (Promega). The HepG2 cells transfected with METAFECTENE showed a much higher activity of β -galactosidase and therefore a much higher transfection efficiency than the cells transfected with the other methods.




Fig. 1: Transfection of HepG2 cells with different methods

In the next step the optimisation protocol for 12-well dishes was performed to transfect the HepG2 cells with METAFECTENE. For this purpose the corresponding amounts of DNA and volumes of METAFECTENE were diluted in serum- and antibiotics-free RPMI medium, mixed, incubated for 20 minutes at room temperature and then added to 1.5×10^5 HepG2 cells in 1 ml serum- and antibiotics-free RPMI medium. After a cultivation period of 24 hours in the incubator the transfection efficiency was determined with the β -Gal assay by the activity of the β -galactosidase.




Fig. 2: Optimisation of the transfection of HepG2 cells with METAFACTENE

The best transfection efficiency was achieved with 1.5 μg pSV β-Gal plasmid and 4 μl METAFACTENE. In the transfections with 0.5 and 1.0 μg plasmid DNA best results were observed with 2 respectively 4 μl METAFACTENE indicating an optimal relation of METAFACTENE to DNA of 4. For this reason a further improvement of the efficiency is probable achievable by the use of 6 μl METAFACTENE in combination with 1,5 μg DNA (METAFACTENE : DNA = 4).