

DNA-transfection of COS 7 and MEF cells using Biontex K2 Transfection System in comparison with Invitrogen Lipofectamine Plus Reagent.

Marzia Ognibene, Molecular Biology Laboratory, G. Gaslini Institute, Largo Gaslini 5, 16147 Genova, Italy

Materials and Methods

Cell cultures

COS 7 and MEF cells were cultured on 100-mm tissue culture dishes (Falcon), in 10 ml of Dulbecco's Modified Eagle's Medium (DMEM) supplemented with L-glutamine, Penicillin+Streptomycin and 10% fetal bovine serum (FBS).

Cells were grown to 70% confluency in transfection test n°1 or to 90% confluency in transfection test n°2. They were transiently transfected with 6 µg of pcef|GST-human-oncoDbl construct (8100 bp).

Cell transfection

Cells were transfected using Lipofectamine Plus (Invitrogen) as described by the manufacturer, or using Biontex K2 Transfection System in parallel.

In the transfection test n°1 we used 300 µl Multiplier and 70 µl of Transfection Reagent, while in the transfection test n°2 we used 200 µl Multiplier and 30 µl of Transfection Reagent. Both transfection tests were performed as follows:

2 hours incubation of cells with Multiplier

Preparation of Solution A with 750 µl serum-free DMEM and 6 µg DNA

Preparation of Solution B with 750 µl serum-free DMEM and Transfection Reagent

Mix of solution A into solution B and incubation of 20-30 minutes

Addition of mixed solution to cells and 24 hours incubation


Replacement with fresh complete growth medium.

Western blot analysis

24 hours after the medium replacement, cells were lysed on ice with a specific lysis buffer and 100 µg of each lysate were subjected to 6% SDS-PAGE electrophoresis, transferred to PVDF membrane (Millipore) and probed with the polyclonal anti-human Dbl antibody (Santa Cruz Biotechnology).

Results and conclusions

Cells from transfection test n°1 showed a certain suffering, especially COS 7, while cells from transfection test n°2 were comparable to cells from Lipofectamine Plus transfection. The Dbl protein expression was surely stronger with K2 Transfection System than with Lipofectamine Plus Reagent, and better with a higher cell confluency combined to a lower amount of K2 reagents (see the figure).



Lane 1: untransfected COS 7 cells

Lane 2: COS 7 cells transfected with Lipofectamine Plus Reagent

Lane 3: COS 7 cells transfected with K2 . Transfection test n°2

Lane 4: untransfected MEF cells

Lane 5: MEF cells transfected with Lipofectamine Plus Reagent

Lane 6: MEF cells transfected with K2 . Transfection test n°1

Lane 7: MEF cells transfected with K2 . Transfection test n°2