

K2 Transfection Technical Note

Clarissa Rocha, Department of Microbiology, University of Sao Paulo, Sao Paulo, Brazil.

Materials: Plasmid: pHIV-dTomato;
DMEM supplement with 10% fetal calf serum and 1% antibiotics;
Cell line: U87MG (glioma).

Protocol:

1×10^4 or 2×10^4 U87MG were plated in 96 well plate with 100 μ L completed medium and incubated in 5% CO₂ at 37°C for 24 hrs. Subsequently, it was performed the transfections as following:

- 1) K2 multiplier and K2 transfection reagent were placed at room temperature and gently mixed;
- 2) 1 μ L of K2 Multiplier was added to each well and incubated for 2 hrs;
- 3) Solution A: 400 ng DNA in 20 μ L
Solution B: **(1:3)** 0,6 μ L K2 Transfection Reagent in 10 μ L serum-free medium
Solution B: **(1:4)** 0,8 μ L K2 Transfection Reagent in 10 μ L serum-free medium;
- 4) 10 μ L of Solution A was added to Solution B and incubated at RT for 15 min.
- 5) 20 μ L of DNA-lipid complex was added to each well. The solution was gently mixed and the cells were incubated overnight at 37°C in CO₂ incubator.
- 6) In the next day, the medium was replaced and transfection efficiency was analyzed in fluorescent microscopy;

RESULTS

In terms of cell viability all tested conditions were not toxic to cells.

In terms of transfection efficiency the best condition was = **1:4** and **2×10^4** cells

Transfection efficiency: U87MG = about **40-50%**

