

- **DNA-transfection of primary human mesenchymal stem cells (hMSCs) using “Biontex K2<sup>®</sup> Transfection System”.**

Shan-hui Hsu, Institute of Polymer Science and Engineering, National Taiwan University, No. 1, Sec. 4 Roosevelt Road, Taipei 10617, Taiwan.

## **Materials and Methods**

### **Cell culture**

Human mesenchymal stem cells (hMSCs) were cultured in  $\alpha$ -MEM supplemented with 10% FBS, 1% penicillin-streptomycin solution, at 37 °C in a humidified atmosphere containing 5% CO<sub>2</sub>. Cells of passage seven were used for the transfection studies.

### **Cell transfection**

5×10<sup>4</sup> human hMSCs were plated in a single well of a 24-well dish in 1 ml medium. On the day of transfection, the cell confluence was 70-90%. Two hours before DNA transfection, 15  $\mu$ l K2<sup>®</sup> Multiplier was dripped to each well. The following solutions were prepared in a polypropylene vessel.

Solution A: 37.5  $\mu$ l serum-free medium with 1  $\mu$ g DNA plasmid.

Solution B: 37.5  $\mu$ l serum-free medium with 2  $\mu$ l K2<sup>®</sup> Transfection Reagent.

Solution A was added to Solution B. The mixture was incubated at room temperature for 20 min. The mixture was then immediately added to the cells. Transfections were conducted at 37 °C and 5% CO<sub>2</sub> for 24 hours. Cell viability and transfection efficiency were estimated.

Dish size	$\alpha$ -MEM	K2 <sup>®</sup> Multiplier	K2 <sup>®</sup> Transfection Reagent	Serum-free medium for DNA/K2 reagent	DNA plasmid
24 well	1 ml	15 $\mu$ l	2 $\mu$ l	37.5 $\mu$ l	1 $\mu$ g

## **Results**


### **Viability of hMSCs**

To determine the cell survival rate, cell pellets were resuspended in PBS containing 1  $\mu$ l of propidium iodide (PI) and subjected to flow cytometric analysis. The viability of hMSCs transfected by K2 transfection reagent was higher than that by a dendrimer-based commercial transfection reagent (Figure 1A).

### **Transfection efficiency of hMSCs**

Real-time RT-PCR was used to analyze the expression of exogenous (transfected)

gene. The mRNA expression of the transfected gene by K2 was about 1.5-fold higher than that by a dendrimer-based commercial transfection reagent (Figure 1B).



### Conclusions

The data shows that K2<sup>®</sup> Transfection System has low cytotoxic effects compared to another dendrimer-based commercial transfection reagent. Besides, the transfection efficiency for hMSCs evaluated by mRNA expression is also higher by K2<sup>®</sup> Transfection System.