

## **DNA transfection of renal epithelial cells, human primary cells and human cell lines, using “Biontex K2<sup>®</sup>Transfection System”.**

Margot Rehm, Margarete Goppelt-Struebe, Nephrologische Forschungslaboratorien,  
Medizinische Klinik 4 - Nephrologie und Hypertensiologie, FAU Erlangen

### **Materials**

K2<sup>®</sup> Multiplier (Biontex Laboratories)

K2<sup>®</sup> Transfection reagent (Biontex Laboratories)

Opti-MEM<sup>®</sup> Reduced Serum Medium (Life Technologies)



### **Cell culture and transfection**

Primary human tubular epithelial cells were isolated and cultured as described previously [1]. The human proximal tubular cell line HKC-8 was kindly provided by S. Racusen [2] and cultured as described [3]. Primary cells (100 000 cells) or HKC8 cells (80 000 cells) were seeded in a 24 well plate on collagen IV-coated glass cover slips in serum containing medium to allow adherence. The next day, medium was changed to 0 % serum, corresponding to the normal growth medium of epithelial cells. Concomitantly, K2<sup>®</sup> Multiplier (5  $\mu$ l / 0.5 ml) was added to the wells. After 2 h, the complex of DNA and K2<sup>®</sup> transfection reagent was added, which had been prepared 20 min before: in a 96 well plate, DNA (0.5  $\mu$ g) was diluted in 15  $\mu$ l Opti-MEM<sup>®</sup> Reduced Serum Medium (Life Technologies); in another well, 15  $\mu$ l of Opti-MEM<sup>®</sup> were mixed with 1.2  $\mu$ l of K2<sup>®</sup> transfection reagent by pipetting once up and down. Subsequently, the diluted DNA was added and mixed by pipetting once up and down. After 15 – 20 min incubation at room temperature the mix was added to the cells (30  $\mu$ l / 500  $\mu$ l). 24 h later, cells were fixed by PFA and proteins of interest were detected by immunocytochemistry using standard techniques.

### References:

1. Kroening S, Neubauer E, Wullich B, Aten J, Goppelt-Struebe M (2010) Characterization of connective tissue growth factor expression in primary cultures of human tubular epithelial cells: modulation by hypoxia. *Am J Physiol Renal Physiol* 298: F796-806.
2. Racusen LC, Monteil C, Sgrignoli A, Lucskay M, Marouillat S, et al. (1997) Cell lines with extended in vitro growth potential from human renal proximal tubule: characterization, response to inducers, and comparison with established cell lines. *J Lab Clin Med* 129: 318-329.
3. Kroening S, Neubauer E, Wessel J, Wiesener M, Goppelt-Struebe M (2009) Hypoxia interferes with connective tissue growth factor (CTGF) gene expression in human proximal tubular cell lines. *Nephrol Dial Transplant* 24: 3319-3325.

## Results

	
<p>Primary tubular epithelial cells were transfected with the transcriptional activator MKL (nuclear localization, red fluorescence); F-actin fibers were visualized with rhodamine-phalloidin (green). Arrows indicate cells with lower expression levels of the transgene.</p>	<p>HKC-8 cells were transfected with dominant negative Rac1 (flag tag detected in green). Endogenous MKL was detected by specific antibodies (pseudo color: red).</p>

Over time, we had tested several transfection reagents most of which were not suitable for transfection of primary cells with transfections efficiencies of 10 – 20 cells per plate. Also with K2<sup>®</sup> transfection reagent transfection efficiencies of primary epithelial cells were not comparable with easy to transfect cell lines, but reached 5-10 %. Viability of the cells was barely affected after 24 h, whereas cells looked less healthy after 48 h. Therefore, we restricted our experiments to 24 h whenever possible.

HKC-8 cells are less difficult to transfect, but also in the lower range of about 10 %. The cells tolerated the transfection very well and were used after 24 without any signs of loss of viability or stress. Transfected HKC-8 cells were also used in reporter gene assays.

Transfection efficiency was very reliable based on the enzyme activity of the control plasmid beta-galactosidase:  $0.504 \pm 0.007$  OD, means  $\pm$  SD of 7 wells transfected in parallel.